



LOYAL[®]

...the soul of surgery...

- Scissors
- Artery Forceps
- Needle Holders
- Micro Scissors
- Micro Needle Holders
- Micro Forceps
- Clamps



CE | ISO | GMP

Welcome

Profile

Quality product is the foundation of **LOYAL**.

LOYAL have been overwhelmed with the response of maxillofacial, plastic & neuro surgeons for our products serving to almost all leading surgeons, hospital, dental & medical college.

Surgeons observed that **LOYAL Products** are match to international standard, such precision work is possible on advance mechanism on CNC, VNC machinery which can maintain zero-zero accuracy.

Our manufacturing setup is capable bulk & huge export order. We have good network in all major cities of India, all raw material is biocompatible medical grade.

LOYAL believes in long term good relation.

Medical device quality management system certified organization.

Thank You!

LOYAL[®]

/// Scissors ///

Iris Scissor - Fine



Super Cut
1501.1
Straight & Curved
4" & 5"



T.C. Tip
1501.1
Straight & Curved
4" & 5"



S.S.
1501.1
Straight & Curved
4" & 5"

Metzenbaum Scissor



Super Cut
1502.1
Straight & Curved
4", 5", 6", 7", 8"



T.C. Tip
1502.2
Straight & Curved
4", 5", 6", 7", 8"



S.S.
1502.3
Straight & Curved
4", 5", 6", 7", 8"

Tenotomy Scissor - Serrated



Super Cut
1503.1
Straight & Curved
4", 5", 6"



T.C. Tip
1503.2
Straight & Curved
4", 5", 6"



S.S.
1503.3
Straight & Curved
4", 5", 6"

/// Scissors ///

Kilner Dissecting Scissor



T.C. Tip
1504.1 - Curved
5"



S.S.
1504.2 - Curved
5"

Mayo Scissor



Super Cut
1505.1
Straight & Curved
5", 6", 7", 8"



T.C. Tip
1505.2
Straight & Curved
5", 6", 7", 8"



S.S.
1505.3
Straight & Curved
5", 6", 7", 8"

Potts - Smith Scissor



T.C. Tip
1506.1 - 45° & 60°
7"



S.S.
1506.2 - 45° & 60°
7"

/// Scissors/ Needle Holders ///

Heath's Stitch Scissor



T.C. Tip
1507.1
Curved to Side
5"



T.C. Tip
1507.2 - Curved
4" & 6"



S.S.
1507.3 - Curved
4" & 6"

Gillies Needle Holder cum Scissor



T.C. Tip
1508.1 - Curved
5.5"



S.S.
1508.2 - Curved
5.5"

Needle Holder with Tungsten Carbide



Mayo (Debakey)
1509.1 - Straight
4", 5", 6", 7"



Bozeman
1509.2 - Angular
6"



Ryder
1509.3 - Straight
4", 5", 6", 7"

/// Scissors/ Needle Holders/ Forceps ///

Needle Holder S.S. - Mayo



Mayo (Debakey)
1510.1 - Straight
4", 5", 6", 7"



Bozeman
1510.2 - Angular
6"



Ryder
1510.3 - Straight
4", 5", 6", 7"

Adson Forceps - Fine



Plain
1511.1
5", 6", 7"



Tooth
1511.2
5", 6", 7"



Plain T.C. Tip
1512.1
5", 6", 7"



Tooth T.C. Tip
1512.2
5", 6", 7"

Dressing Forceps



Plain
1513.1
5", 6", 7"



Tooth
1513.2
5", 6", 7"



Debakey
1514
5", 6", 7"

/// Artery/ Needle Holders/ Forceps ///

Bayonet Forceps



Plain
1515.1
6", 7"

Tooth
1515.2
6", 7"

Mosquito Artery Forceps



Straight
1516.1
4", 5"



Curved
1516.2
4", 5"

Crile Artery Forceps



Straight
1517.1
6", 7", 8"



Curved
1517.2
6", 7", 8"

Kocher Artery Forceps



Straight
1518.1
6", 8", 10"



Curved
1518.2
6", 8", 10"

/// Artery/ Forceps/ Towel Clips ///

Mixer Artery Forceps



1519
6", 7", 8"

Sponge Holding Forceps



1520
6", 8", 10"

Allis Tissue Forceps



1521
6", 8", 10"

Babcock Tissue Forceps



1522
6", 8", 10"

Towel Clips



Backhaus Forcep
1523.1



'X' Action Forcep
1523.2

/// B.P. Handles/ Forceps ///

B.P. Handles



Round Handle
1524.1
No. 3, No. 4



Scalpel Handle D
1524.2
No. 7



Scalpel Handle D/E
1524.3
No. 3, No. 4



Scalpel Handle
1524.4
No. 3, No. 4

Skin Staple Removal



1525

Cheatle Forceps



1526

/// Micro Scissors ///

Micro Scissors S.S.



Straight
1550.1



Curved
1550.2

/// Micro Scissors/ Needle Holders/ Forceps ///

Micro Needle Holder S.S.



Straight & Curved
1551

Micro Scissors Titanium



Straight & Curved
1552

Micro Needle Holder Titanium



Straight & Curved
1553

Micro Bayonet Scissor




1554
6"

Micro Forceps S.S.



Plain Tooth
1555
5", 6"



Plain Tooth
1556
5", 6"

Micro Forceps / Clamps

Micro Forceps Titanium



Plain Tooth
1557
5", 6"

Plain Tooth
1558
5", 6"

Micro Forceps Ring Tip



1559
5", 6"

Bull dog Clamps



Straight & Curved
1560
38mm/1", 51mm/2", 57mm/2.1/4", 63mm/2.1/2", 75mm/3

Jeweller Forceps



Straight & Curved
1561
38mm/1", 51mm/2", 57mm/2.1/4", 63mm/2.1/2", 75mm/3


Vessel Dialater



Straight & Curved
1562
38mm/1", 51mm/2", 57mm/2.1/4", 63mm/2.1/2", 75mm/3

/// Micro Forceps ///

Tweezer Forceps



1563

/// Kits ///

GENERAL SURGERY KIT

CONTENTS	QTY
Iris Scissors - Fine Straight	01 Pcs
Iris Scissors -Curved	01 Pcs
Adson Dressing Forceps Plain	01 Pcs
Adson Dressing Forceps Tooth	01 Pcs
Allis Tissue Forcep 6"	02 Pcs
Backhaus Towel Clips	02 Pcs
Crile Artery Forcep Straight 6 ½"	01 Pcs
Crile Artery Forcep Curved 6 ½"	01 Pcs
Sponge holding8"	01Pcs
Sponge holding10"	01 Pcs
Mosquito Artery Forcep Straight 5"	02 Pcs
Mosquito Artery Forcep Curved 5"	02 Pcs
Mayo Dissecting Scissors Straight	01 Pcs
Metzenbaum Scissors Curved 6"	01 Pcs
Suction Cannula - MastoidNo.3	01 Pcs
B P Handle - RoundNo.3	01 Pcs
Health's Stitch Scissors	01 Pcs
Needle Holder With T.C.Tip 5"	01 Pcs
Needle Holder With T.C.Tip 6"	01 Pcs
Kilner Dissecting Scissors5"	01 Pcs
Tenotomy Scissors Curved 5"	01 Pcs
Potts - Smith Scissors	01 Pcs
Towel Clip X Action	01 Pcs
Aluminium Box	01 Pcs

/// Kits ///

INSTRUMENT SET FOR MICRO SURGERY KIT

CONTENTS	QTY
Needle Holders Castroviejo - 6" Straight	01 Pcs
Needle Holders Castroviejo - 6" Curved	01 Pcs
Micro Forceps 6" Plain	01 Pcs
Micro Forceps 6" Platform	01 Pcs
Scissors - Castroviejo Straight	01 Pcs
Scissors - Castroviejo Curved	01 Pcs
Jeweller's Forceps	01 Pcs
General Instruments Box with Silicone Mate - Small	01 Pcs

MINOR ORAL SURGERY KIT

CONTENTS	QTY
Adson - Forceps Fine - PLAIN 7"	01 Pcs
Adson - Forceps Fine - TOOTH 7"	01 Pcs
Allis Tissue Forceps 6"	02 Pcs
Crile Artery Forceps -St.6"	01 Pcs
Crile Artery Forceps -Cd .6"	01 Pcs
Sponge Holding Forceps - 8"	01 Pcs
Sponge Holding Forceps - 10"	01 Pcs
Mosquito Artery Forceps - St.4.5"	02 Pcs
Mosquito Artery Forceps - Cd.4.5"	02 Pcs
Mayo Scissors - Cd.6"	01 Pcs
Metzenbaum Scissors -Cd.7"	01 Pcs
Suction Cannula - Mastoid No.3	01 Pcs
B P Handle - Round Handle	01 Pcs
Heath's Stitch Scissors	01 Pcs
Needle Holder with Tungsten Carbide	02 Pcs
Tenotomy Scissors -Cd.5"	01 Pcs
Dressing Forceps -Plain	01 Pcs
Poswillo Malar Hook	01 Pcs
Aluminium Instruments Box	01 Pcs

Certificate

Certificate



Royal Stancert Inc.
★★★★★★★★★★★★

This is to Certify that

LOYAL

FF 33, Harmony Commercial Mall, New Link Road, Bhagat Singh Nagar - 1, Goregaon (West), Mumbai - 400104 (Maharashtra), India.

has been assessed by RSI and found to comply with the requirements of

ISO 13485:2016
Medical Device - Quality Management System

For the following Scope:
Manufacturing and Trading of Maxillofacial, Plastic, Neuro, Orthopaedic and General Surgery Implants & Instruments.

Certificate No.: M-RS-1803-3174
Date of Issue: 30-03-2018
Date of Expiry: 29-03-2021
1st Surveillance due: 29-02-2019
1ind Surveillance due: 29-02-2020




Director (Certification)




PAC-GEAC-1506-299

International Operations: Eschersheimer Landstr. 42, 60322 Frankfurt am Main, Frankfurt, Germany.
Email: info@royalstancert.net

This certificate remains the property of Royal Stancert Inc. and must be returned whenever demanded. The validity of this certificate can be verified at <http://www.royalstancert.net>. Royal Stancert Inc. is an independent system, product and personal assessment body. RSI is accredited by Global Euro Accreditation Centre, Georgia. (GCN - 654)

Certificate



Royal Stancert Inc.
★★★★★★★★★★★★

This is to Certify that

LOYAL

FF 33, Harmony Commercial Mall, New Link Road, Bhagat Singh Nagar - 1, Goregaon (West), Mumbai - 400104 (Maharashtra), India.

has been assessed by RSI and found to comply with the requirements of

ISO 9001:2015
Quality Management System

For the following Scope:
Manufacturing and Trading of Maxillofacial, Plastic, Neuro, Orthopaedic and General Surgery Implants & Instruments.

Certificate No.: Q-RS-1803-3173
Date of Issue: 30-03-2018
Date of Expiry: 29-03-2021
1st Surveillance due: 29-02-2019
1ind Surveillance due: 29-02-2020




Director (Certification)




PAC-GEAC-1506-299

International Operations: Eschersheimer Landstr. 42, 60322 Frankfurt am Main, Frankfurt, Germany.
Email: info@royalstancert.net

This certificate remains the property of Royal Stancert Inc. and must be returned whenever demanded. The validity of this certificate can be verified at <http://www.royalstancert.net>. Royal Stancert Inc. is an independent system, product and personal assessment body. RSI is accredited by Global Euro Accreditation Centre, Georgia. (GCN - 654)

Certificate



Royal Stancert Inc.
★★★★★★★★★★★★

This is to certify that the technical documentation for the product:

Maxillofacial, Plasic, Neuro and General Surgery Implants and Instruments.
(Product Details as per Annexure 1)

Manufacturing By:
LOYAL

FF 33, Harmony Commercial Mall, New Link Road, Bhagat Singh Nagar - 1, Goregaon (West), Mumbai - 400104 (Maharashtra), India.

Complies with the applicable requirements of **Medical Device Directive 93/42/EEC**

The technical documentation / inspection / test results comply with the requirements of **Medical Device Directive (MDD)** hence manufacturer's declaration of conformity according to above directive is accepted **Loyal** can place the CE marking as per laid down regulations on the Products mention as per Annexure.

Certificate No.: CE-RS-1803-3176
Date of Issue: 30-03-2018
Date of Expiry: 29-03-2021
Surveillance Frequency: Annual




Director (Certification)




PAC-GEAC-1506-299

International Operations: Eschersheimer Landstr. 42, 60322 Frankfurt am Main, Frankfurt, Germany.
Email: info@royalstancert.net

This certificate remains the property of Royal Stancert Inc. and must be returned whenever demanded. The validity of this certificate can be verified at <http://www.royalstancert.net>. Royal Stancert Inc. is an independent system, product and personal assessment body. RSI is accredited by Global Euro Accreditation Centre, Georgia. (GCN - 654)

Certificate



Royal Stancert Inc.
★★★★★★★★★★★★

This is to Certify that

LOYAL

FF 33, Harmony Commercial Mall, New Link Road, Bhagat Singh Nagar - 1, Goregaon (West), Mumbai - 400104 (Maharashtra), India.

has been assessed by RSI and found to comply with the requirements of

GMP
(Guidelines for Good Manufacturing Practice)

For the following Scope:
Manufacturing and Trading of Maxillofacial, Plastic, Neuro, Orthopaedic and General Surgery Implants & Instruments.

Certificate No.: G-RS-1803-3175
Date of Issue: 30-03-2018
Date of Expiry: 29-03-2021
1st Surveillance due: 29-02-2019
1ind Surveillance due: 29-02-2020




Director (Certification)




PAC-GEAC-1506-299

International Operations: Eschersheimer Landstr. 42, 60322 Frankfurt am Main, Frankfurt, Germany.
Email: info@royalstancert.net

This certificate remains the property of Royal Stancert Inc. and must be returned whenever demanded. The validity of this certificate can be verified at <http://www.royalstancert.net>. Royal Stancert Inc. is an independent system, product and personal assessment body. RSI is accredited by Global Euro Accreditation Centre, Georgia. (GCN - 654)

CHECK THE HARDNESS OF WATER BEFORE USE.
DO NOT USE BORING WATER OR WHICH HAS HIGH DEPOSITS OF SALT IN IT

Care and Maintenance

- ALWAYS** Check all instruments before use, Unless otherwise directed, they should be lubricated, cleaned and washed before being put into sterilisation.
- ALWAYS** Ensure instruments are only used for the purpose for which they have been designed.
- ALWAYS** Handle all instruments gently, never overstrain, drop or misuse them.
- ALWAYS** Check all instruments for damage after use.
- ALWAYS** Clean and decontaminate all instruments in cold or warm water as soon as possible after use. Failure to do so may result in the instruments becoming stained and the joints stiff.
- ALWAYS** Give special attention to micro-surgical instruments. Their fine tips can easily be damaged by contact with other instruments or the sides of the case in which they are kept.
- ALWAYS** Lubricate after cleaning with a proprietary water soluble lubricant.
- ALWAYS** Ensure all instruments are thoroughly dried immediately after they are washed. Failure to do so the instruments will get stained.
- ALWAYS** While autoclaving pack all heavy instruments lying on a piece of cloth or towel at the bottom, to avoid all instruments getting in direct contact with each other.
- ALWAYS** Store and sterilise bow handled instruments on a special holder, rack and rachets left open.
- ALWAYS** Check the hardness of the water in the autoclave, too hard water will leave deposits on the instruments.
- ALWAYS** Return all instruments to us for replacement or repair in guarantee period
- ALWAYS** Replace all rubber washers after one year.
-
- NEVER** Leave soiled instruments to dry. If it is impossible to clean them immediately after use, then soak them in water for as short period as possible.
- NEVER** Overstrain joints or rachets.
- NEVER** Misuse surgical instruments for any other use.
- NEVER** Use abrasives on instruments as this will spoil the surface finish. This may later cause de-colouration, rusting or pitting.
- NEVER** Leave instruments soaking longer than necessary in chemical sterilising solutions, after they have been used, wash thoroughly in warm water to remove all traces of chemicals.
- NEVER** Store damp instruments, they must be thoroughly dried before storing

Thank you.

LOYAL®

Office:

FF 33, Harmony Commercial Mall,
Bhagat Singh Nagar No. 1, New Link Road,
Goregaon (West), Mumbai - 400104,
Contact: 91-9821225006 / 9322941234.
Email: loyalsurgical@gmail.com
info@loyalsurgical.com



LOYAL[®]
...the soul of surgery...

Suggestion & Complain No.
91-9702484003.

Factory:

Somnath Ortho Care Pvt. Ltd, GIDC,
Valsad, Umbergaon, Gujarat - 396171

www.loyalsurgical.com